

Insert P/N 4425879 Rev. C  
Printed in USA

For Research Use Only.  
Not for use in diagnostic procedures.

Part Number	Description	Volume	Storage Conditions
4415204	VetMAX™ BRSV Reagents	110 µL	-20 °C. <b>Do not store in a frost-free freezer.</b>

**Note:** Read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

## Product description

The VetMAX™ BRSV Reagents is a primer and probe mix designed to amplify bovine respiratory syncytial virus (BRSV) RNA and Xeno™ RNA Control in multiplex real-time RT-PCR. The primer and probe designs for BRSV RNA amplification are published and publicly available.

Use 1 µL of the VetMAX BRSV Reagents per 25-µL reaction.

## TaqMan® Probe Reporter/Quencher

Target	Reporter	Quencher
BRSV RNA	FAM™ dye	Eclipse® Q
Xeno™ RNA Control	VIC® dye	Eclipse® Q

## Safety information

### Obtaining MSDSs

To obtain Material Safety Data Sheets (MSDSs) for any chemical product supplied by Applied Biosystems or Ambion:

- At [www.appliedbiosystems.com](http://www.appliedbiosystems.com), select **Support** ▶ **MSDS**. Search by chemical name, product name, product part number, or MSDS part number. Right-click to print or download the MSDS of interest.
- At [www.ambion.com](http://www.ambion.com), go to the web catalog page for the product of interest. Select **MSDS**, then right-click to print or download.
- E-mail (MSDS\_Inquiry\_CCRM@appliedbiosystems.com), telephone (650-554-2756; USA), or fax (650-554-2252; USA) your request, specifying the catalog or part number(s) and the name of the product(s). The associated MSDSs will be e-mailed unless you request fax or postal delivery. Requests for postal delivery require 1 to 2 weeks for processing.

**Note:** For the MSDSs of chemicals not distributed by Applied Biosystems or Ambion, contact the chemical manufacturer.

---

For Research Use Only. Not for use in diagnostic procedures.

Information in this document is subject to change without notice.

APPLIED BIOSYSTEMS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL APPLIED BIOSYSTEMS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT APPLIED BIOSYSTEMS IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

**Warranty and Liability:** Applied Biosystems is committed to delivering superior product quality and performance, supported by industry-leading global service and technical support teams. Warranty information for the accompanying consumable product is available at [www.ambion.com/info/warranty](http://www.ambion.com/info/warranty) in "Limited Warranty for Consumables," which is subject to the exclusions, conditions, exceptions, and limitations set forth under the caption "EXCLUSIONS, CONDITIONS, EXCEPTIONS, AND LIMITATIONS" in the full warranty statement. Please contact Applied Biosystems if you have any questions about our warranties or would like information about post-warranty support.

Limited Use Label License: Veterinary Diagnostics

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform veterinary diagnostic services including reporting the results of purchaser's activities for a fee and for the purchaser's own internal research. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008. No right under any other patent claim (such as claims to methods, apparatus or reagents) is conveyed expressly, by implication, or by estoppel.

**Trademarks:** Trademarks of Life Technologies Corporation and its affiliated companies: AB Logo™, Applied Biosystems®, FAM™, VetMAX™, VIC®, and Xenon™. TaqMan is a registered trademark of Roche Molecular Systems, Inc. All other trademarks are the sole property of their respective owners.

© Copyright 2009, 2011 Life Technologies Corporation. All rights reserved.

Part Number 4425879 Rev. C 06/2011

---



**Headquarters**

850 Lincoln Centre Drive | Foster City, CA 94404 USA  
Phone 650.638.5800 | Toll Free 800.345.5224  
[www.appliedbiosystems.com](http://www.appliedbiosystems.com)

**International Sales**

For our office locations please call the division headquarters or refer to our Web site at [www.appliedbiosystems.com/about/offices.cfm](http://www.appliedbiosystems.com/about/offices.cfm)