

Part Number	Description	Volume	Storage Conditions
4415192	VetMAX™ BRSV Controls	80 µL	-20 °C. <i>Do not store in a frost-free freezer.</i>

Note: Read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product description

The VetMAX™ BRSV Controls is a mixture of bovine respiratory syncytial virus (BRSV) control RNA (1000 copies/µL) and Xeno™ RNA Control (1000 copies/µL) that serves as a positive control for amplification of BRSV RNA and Xeno RNA Control in multiplex real-time RT-PCR.

Use 8 µL of the VetMAX BRSV Controls (equivalent to 8000 copies each of BRSV control RNA and Xeno RNA Control) per 25-µL reaction.

Handling instructions

Use typical precautions against RNase contamination when handling the VetMAX BRSV Controls. For example, wear clean gloves and use nuclease-free barrier pipet tips. Improper handling of the VetMAX BRSV Controls can result in RNA degradation and failure of the positive control RT-PCR.

Safety information

Obtaining MSDSs

To obtain Material Safety Data Sheets (MSDSs) for any chemical product supplied by Applied Biosystems or Ambion:

- At www.appliedbiosystems.com, select **Support** ▶ **MSDS**. Search by chemical name, product name, product part number, or MSDS part number. Right-click to print or download the MSDS of interest.
- At www.ambion.com, go to the web catalog page for the product of interest. Select **MSDS**, then right-click to print or download.
- E-mail (MSDS_Inquiry_CCRM@appliedbiosystems.com), telephone (650-554-2756; USA), or fax (650-554-2252; USA) your request, specifying the catalog or part number(s) and the name of the product(s). The associated MSDSs will be e-mailed unless you request fax or postal delivery. Requests for postal delivery require 1 to 2 weeks for processing.

Note: For the MSDSs of chemicals not distributed by Applied Biosystems or Ambion, contact the chemical manufacturer.

For Research Use Only. Not for use in diagnostic procedures.

Information in this document is subject to change without notice.

APPLIED BIOSYSTEMS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL APPLIED BIOSYSTEMS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT APPLIED BIOSYSTEMS IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Warranty and Liability: Applied Biosystems is committed to delivering superior product quality and performance, supported by industry-leading global service and technical support teams. Warranty information for the accompanying consumable product is available at www.ambion.com/info/warranty in "Limited Warranty for Consumables," which is subject to the exclusions, conditions, exceptions, and limitations set forth under the caption "EXCLUSIONS, CONDITIONS, EXCEPTIONS, AND LIMITATIONS" in the full warranty statement. Please contact Applied Biosystems if you have any questions about our warranties or would like information about post-warranty support.

Limited Use Label License: Veterinary Diagnostics

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform veterinary diagnostic services including reporting the results of purchaser's activities for a fee and for the purchaser's own internal research. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008. No right under any other patent claim (such as claims to methods, apparatus or reagents) is conveyed expressly, by implication, or by estoppel.

Trademarks: The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

© Copyright 2010, 2011 Life Technologies Corporation. All rights reserved.

Part Number 4425884 Rev. C 06/2011
