

Insert P/N 4425901 Rev.C  
Printed in USA

For Research Use Only.  
Not for use in diagnostic procedures.

Part Number	Description	Volume	Storage Conditions
4415221	VetMAX™ <i>T. foetus</i> Reagents	110 µL	-20 °C. <b>Do not store in a frost-free freezer.</b>

**Note:** Read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

## Product description

The VetMAX™ *T. foetus* Reagents is a primer and probe mix designed to amplify *Tritrichomonas foetus* DNA and Xeno™ DNA Control in multiplex real-time PCR. The primer and probe designs for *T. foetus* DNA amplification are published and publicly available.

Use 1 µL of the VetMAX *T. foetus* Reagents per 25-µL reaction.

## TaqMan® Probe Reporter/Quencher

Target	Reporter	Quencher
<i>T. foetus</i> DNA	FAM™ dye	MGB
Xeno™ DNA Control	VIC® dye	Eclipse® Q

## Safety information

### Obtaining MSDSs

To obtain Material Safety Data Sheets (MSDSs) for any chemical product supplied by Applied Biosystems or Ambion:

- At [www.appliedbiosystems.com](http://www.appliedbiosystems.com), select **Support** ▶ **MSDS**. Search by chemical name, product name, product part number, or MSDS part number. Right-click to print or download the MSDS of interest.
- At [www.ambion.com](http://www.ambion.com), go to the web catalog page for the product of interest. Select **MSDS**, then right-click to print or download.
- E-mail (MSDS\_Inquiry\_CCRM@appliedbiosystems.com), telephone (650-554-2756; USA), or fax (650-554-2252; USA) your request, specifying the catalog or part number(s) and the name of the product(s). The associated MSDSs will be e-mailed unless you request fax or postal delivery. Requests for postal delivery require 1 to 2 weeks for processing.

**Note:** For the MSDSs of chemicals not distributed by Applied Biosystems or Ambion, contact the chemical manufacturer.

For Research Use Only. Not for use in diagnostic procedures.

Information in this document is subject to change without notice.

APPLIED BIOSYSTEMS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL APPLIED BIOSYSTEMS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT APPLIED BIOSYSTEMS IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

**Warranty and Liability:** Applied Biosystems is committed to delivering superior product quality and performance, supported by industry-leading global service and technical support teams. Warranty information for the accompanying consumable product is available at [www.ambion.com/info/warranty](http://www.ambion.com/info/warranty) in "Limited Warranty for Consumables," which is subject to the exclusions, conditions, exceptions, and limitations set forth under the caption "EXCLUSIONS, CONDITIONS, EXCEPTIONS, AND LIMITATIONS" in the full warranty statement. Please contact Applied Biosystems if you have any questions about our warranties or would like information about post-warranty support.

Limited Use Label License: Veterinary Diagnostics

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform veterinary diagnostic services including reporting the results of purchaser's activities for a fee and for the purchaser's own internal research. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008. No right under any other patent claim (such as claims to methods, apparatus or reagents) is conveyed expressly, by implication, or by estoppel.

**Notice To Purchaser: Limited Licenses**

Patents covering oligonucleotide conjugates of Minor Groove Binder ("MGB") are owned by Nanogen and licensed exclusively to Applied Biosystems for use in assays whereby the detection is mediated by a 5' nuclease activity of a polymerase enzyme (the "5' Nuclease Process"). TaqMan® MGB Probes are sold by Applied Biosystems solely for the purchaser's own use in performing 5' Nuclease assays and not for repackaging or resale in any form including resale by purchasers who are licensed to make and sell kits for use in the 5' Nuclease Process. The use of TaqMan® MGB Probes in diagnostics requires separate licenses from Nanogen and from Roche. Use of the 5' Nuclease Process in performing commercial services of any kind (other than diagnostics), including reporting the results of purchaser's activities for a fee requires a separate license from Applied Biosystems. Further information may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

**Trademarks:** Trademarks of Life Technologies Corporation and its affiliated companies: AB Logo™, Applied Biosystems®, FAM™, VetMAX™, VIC®, and Xeno™. TaqMan is a registered trademark of Roche Molecular Systems, Inc. All other trademarks are the sole property of their respective owners.

© Copyright 2010, 2011 Life Technologies Corporation. All rights reserved.

Part Number 4425901 Rev. C 06/2011

---