

RapidMAX™ Software

Create and analyze a RapidMAX™ file

4469875 Rev. A Revision Date September 2011

This quick reference card covers:

■ Log in and install assays	1
■ Import data and set up for analysis	1
■ Set up data for analysis	1
■ Analyze the data	2
■ View analysis results	2
■ Export results	3

Log in and install assays

1. Double-click the RapidMAX™ Software software icon on your desktop.
2. Log in with your user name and password.

Import data and set up for analysis

1. In the Home screen, click **Import File(s) to Analyze** to import one or more files.
2. Select files from the appropriate folder.
3. Assign plate properties by following the Plate Setup Wizard.
Optional: If you want to enter information on the reagents used or enter a description of the file, click **Properties** in the Setup stage of the workflow.

Set up data for analysis

Note: You can use the Plate Grid or the Well Table to set up data.

- 1 Assign assays
 - a. Select wells:

Note: In the Plate Grid only, use any of the following methods to select wells:

 - Click and drag your mouse across wells.
 - Ctrl-click on specific wells you want to select.
 - Click the top left corner of the Plate Grid (above row A, left of column 1) to select all wells.
 - Click the heading of a row or column.
 - b. Right-click, then assign assays.

-
- 2** Assign tasks (controls or unknowns)
- In the Plate Grid, select positive control wells.
 - Right-click, then select **Assign tasks ▶ POSITIVE CONTROL**.
 - Select no-template control (NTC) wells, right-click, then select **Assign tasks ▶ NTC**.
 - (Optional) Select NEC wells, right-click, then select **Assign tasks ▶ NEGATIVE EXTRACTION CONTROL**.

Note: Refer to *RapidMAX™ Software User Guide* (PN 4469874A) for information about importing sample IDs from Microsoft® Excel®.

-
- 3** (Optional) Assign sample ID
- Click on drop-down menu directly above well A1.
 - Select **Sample**.
 - Select wells, then enter the sample ID in the blank data entry box.
-

Analyze the data

-
- 1** Analyze the RapidMAX™ file
- In the taskbar, click  **Analyze** to generate analysis results.

Note: If the Analyze button is faded (not activated), one or more of the following attributes have not been specified in the plate:

- Assays for each well
- Location of positive controls
- Location of NTCs

-
- 2** Save the results
- Click **Save** to save the results, or click **Save as** to save the results as a new file.

Note: You can save results anytime.

View analysis results

-
- 1** View plate results
- In the Workflow pane, click **Results ▶ Plate Results** to open the Plate Results screen and view well details on a plate-by-plate basis.
 - In the Plate Grid or Well Table, right-click a well.
 - Select an action:
 - **Omit well:** Remove well or sample information from the result report
 - **Change call:** An authorized user can change the call on user-defined assays.
 - In the Plate Grid, click **Export** to image to save an image of the analysis results.

-
- 2** View amplification plots
- To view amplification plots:
- In the Plate Grid tab, double-click a well or wells to open the QC Plots window.
 - Click **Results ▶ Amplification Results**.
-

-
- 3** View plate results in the Well Table and check warning flags
- a. Click the **Well Table** tab.
 - b. Click **View**.
 - c. Select or deselect columns to show or hide columns.
 - d. To sort results, click a column header.
Note: To sort data in multiple columns, press the **Ctrl** key and click column headings.
 - e. To see amplification plots, double-click a well row.
-

- 4** View sample results
- To view results combined across multiple targets (for example, subtyping), click **Sample Results** in the Workflow pane.
-

- 5** View amplification results
- a. In Workflow-Results, click **Amplification Results**.
 - b. Select wells to view amplification and multicomponent plots.
 - c. To view presorted results, click a tab: All Calls, Inconclusive Calls, Positive Calls, or Negative calls.
-

Export results

-
- 1** Export analysis results in table format
- a. In the Workflow pane, select **Export ▶ Export**.
 - b. Select the content that you want to export.
 - c. Browse for the export destination of the file, enter an export file name, choose a file type (*.csv, *.txt, *.xls, or *.xlsx), then click **Start Export**.
 - d. Locate the exported file, then open the file.
-

- 2** Export results as a PDF
- a. In the Workflow pane, select **Export ▶ Report**.
 - b. In the Report Preview screen, select the type of report to export.
 - c. Click **Settings** to define the setups for the page, headers and footers, and plate results table.
 - d. Click **Preview** to view the report.
 - e. To accept the report, click the E-Sig checkbox.
 - f. Click **Export**.
 - g. Locate the exported report (*.pdf file), then select **Open With ▶ Adobe Acrobat** (or Adobe Reader) to view the report.
-

Not intended for any human therapeutic or diagnostic use.

NOTICE TO PURCHASER: PLEASE REFER TO THE RAPIDMAX™ SOFTWARE USER GUIDE (PN 4469874) FOR LIMITED LABEL LICENSE OR DISCLAIMER INFORMATION.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

Excel and Microsoft are registered trademarks of Microsoft Corporation.

© 2011 Life Technologies Corporation Corporation. All rights reserved.

Headquarters

5791 Van Allen Way | Carlsbad, CA 92008 USA | Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

For support visit www.appliedbiosystems.com/support

www.lifetechnologies.com

