

PLASTIC CONSUMABLES COMPATIBILITY CHART FOR APPLIED BIOSYSTEMS® INSTRUMENTS

Category/Part Description	P/N	End-Point PCR Systems								Real-Time PCR Systems							Genetic Analyzers		
		Veriti® 96-Well 0.1 mL	Veriti® 96-Well 0.2 mL	Veriti® 384-Well	60-Well Veriti®/9700 (60-Well)	2720	6100	9700 (96-Well)	9700 (384-Well)	StepOne™	StepOnePlus™	7000	7300/7500	7500 Fast	7900HT (96-Well)	7900HT (384-Well)	7900HT Fast	310*	3100/3130 3700/3730
96-Well Plates																			
MicroAmp® Fast Optical 96-Well Reaction Plate (0.1 mL) - 10 plates	4346907	x								x				x			x	x	x
MicroAmp® Fast Optical 96-Well Reaction Plate with Barcode (0.1 mL) - 20 plates	4346906	x								x				x			x	x	x
MicroAmp® Fast Optical 96-Well Reaction Plate with Barcode (0.1 mL) - 200 plates	4366932	x								x				x			x	x	x
MicroAmp® Optical 96-Well Reaction Plate - 10 plates	N8010560		x			x	x	x			x	x		x					x
MicroAmp® Optical 96-Well Reaction Plate - 500 plates	4316813		x			x	x	x			x	x		x					x
MicroAmp® Optical 96-Well Reaction Plate with Barcode - 20 plates	4306737		x			x	x	x			x	x		x					x
MicroAmp® Optical 96-Well Reaction Plate with Barcode - 500 plates	4326659		x			x	x	x			x	x		x					x
MicroAmp® Optical 96-Well Reaction Plate with Barcode and Optical Caps - 20 plates	403012		x			x		x			x	x		x					
MicroAmp® Optical 96-Well Reaction Plate with Barcode and Optical Adhesive Films - 100 plates	4314320		x					x			x	x		x					
384-Well Plates																			
MicroAmp® Optical 384-Well Reaction Plate with Barcode - 50 plates	4309849			x					x						x				x
MicroAmp® Optical 384-Well Reaction Plate with Barcode - 500 plates	4326270			x					x						x				x
MicroAmp® Optical 384-Well Reaction Plate with Barcode- 1,000 plates	4343814			x					x						x				x
MicroAmp® Optical 384-Well Reaction Plate - 1,000 plates	4343370			x					x						x				x
48-Well Plates																			
MicroAmp® Fast Optical 48-Well Reaction Plate - 20 plates	4375816	x								x									
Single Tubes																			
MicroAmp® Fast Reaction Tube with Cap (0.1 mL) - 1,000 tubes	4358297	x								x	x								
MicroAmp® Reaction Tube with Cap (0.2 mL) - 1,000 tubes	N8010540		x			x		x											
MicroAmp® Reaction Tube with Cap (0.2 mL) - 10,000 tubes	N8011540		x			x		x											
MicroAmp® Reaction Tube with Cap, Assorted Colors (0.2 mL) - 1,000 tubes	N8010840		x			x		x											
MicroAmp® Reaction Tube with Cap, Autoclaved (0.2 mL) - 1,000 tubes	N8010612		x			x		x											
MicroAmp® Reaction Tube without Cap (0.2 mL) - 2,000 tubes	N8010533		x			x		x											
MicroAmp® Reaction Tube without Cap (0.2 mL) - 10,000 tubes	N8011533		x			x		x											
MicroAmp® Reaction Tube without Cap, Assorted Colors (0.2 mL) - 1,000 tubes	N8010833		x			x		x											
MicroAmp® Optical Tube without Cap (0.2 mL) - 2,000 tubes	N8010933		x			x		x				x	x						
GeneAmp® Thin-Walled Reaction Tube with Flat Cap (0.5 mL) - 1,000 tubes	N8010737				x														
GeneAmp® Thin-Walled Reaction Tube with Domed Cap (0.5 mL) - 2,000 tubes	N8010537				x														
GeneAmp® Thin-Walled Reaction Tube with Domed Cap, Autoclaved (0.5 mL) - 1,000 tubes	N8010611				x														
MicroAmp® 96-Well Reaction Tube/Tray/Retainer Set (0.2 mL) - 20 trays	403083					x		x											
MicroAmp® 96-Well Reaction Tube/Tray/Retainer Set (0.2 mL) - 100 trays	403086					x		x											
Strip Tubes and Caps																			
MicroAmp® 12-Cap Strip - 200 strips	N8010534	x	x			x		x											
MicroAmp® 12-Cap Strip - 1,000 strips	N8011534	x	x			x		x											
MicroAmp® 12-Cap Strip, Assorted Colors - 200 strips	N8010834	x	x			x		x	x										
MicroAmp® Fast 8-Tube Strip (0.1 mL) - 125 strips	4358293	x							x	x			x				x		
MicroAmp® 8-Tube Strip (0.2 mL) - 125 strips	N8010580		x			x		x											x
MicroAmp® 8-Tube Strip, Assorted Colors (0.2 mL) - 120 strips	N8010838		x			x		x	x										x
MicroAmp® Optical 8-Tube Strip (0.2 mL) - 125 strips	4316567		x			x		x	x			x	x		x				
MicroAmp® 8-Cap Strip - 300 strips	N8010535	x	x			x		x											
MicroAmp® 8-Cap Strip - 1,500 strips	N8011535	x	x			x		x											
MicroAmp® 8-Cap Strip, Assorted Colors - 300 strips	N8010835	x	x			x		x											
MicroAmp® Optical 8-Cap Strip - 300 strips	4323032	x	x			x		x	x	x	x	x	x	x		x			
Seals and Covers																			
MicroAmp® Clear Adhesive Film - 100 films	4306311	x	x	x				x	x										
MicroAmp® Optical Adhesive Film - 25 films	4360954	x	x	x				x	x	x	x	x	x	x	x	x			
MicroAmp® Optical Adhesive Film - 100 films	4311971	x	x	x				x	x	x	x	x	x	x	x	x			
MicroAmp® 48-Well Optical Adhesive Film - 100 films	4375323	x							x										
MicroAmp® 48-Well Optical Adhesive Film - 25 films	4375928	x							x										
MicroAmp® Optical Adhesive Film Kit - 1 kit	4313663							x	x			x		x					
MicroAmp® 96-Well Full Plate Cover - 5 covers	N8010550							x											
Reaction Trays																			
MicroAmp® 96-Well Tray/Retainer Set - 10 trays	403081					x		x				x							
MicroAmp® 96-Well Tray - 10 trays	N8010541					x		x											
MicroAmp® 96-Well Tray for VeriFlex™ Blocks - 10 trays	4379983	x	x							x									
MicroAmp® 96-Well Tray/Retainer Set for Veriti® Systems - 10 trays	4381850		x																
MicroAmp® Fast 48-Well Tray - 10 trays	4375282								x										
Compression Pads																			
MicroAmp® Optical Film Compression Pad - 5 pads	4312639										x			x					
MicroAmp® Snap-On Optical Film Compression Pad - 9 pads	4333292													x					

* Requires the use of proper adapters.



www.invitrogen.com

For research use only. Not intended for any animal or human therapeutic or diagnostic use, unless otherwise stated.

NOTICE TO PURCHASER: PLEASE REFER TO THE USER'S GUIDE OR PACKAGE INSERT OF THE PRODUCTS NAMED HEREIN FOR LIMITED LABEL LICENSE OR DISCLAIMER INFORMATION.

© 2010 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. CO11983 0310