

PLASTIC CONSUMABLES COMPATIBILITY CHART FOR APPLIED BIOSYSTEMS® INSTRUMENTS

Category/Part Description	P/N	End-Point PCR Systems								Real-Time PCR Systems							Genetic Analyzers		
		Veriti® 96-Well 0.1 mL	Veriti® 96-Well 0.2 mL	Veriti® 384-Well	60-Well Veriti®/ 9700 (60-Well)	2720	6100	9700 (96-Well)	9700 (384-Well)	StepOne™	StepOnePlus™	7000	7300/7500	7500 Fast	ViiA™ 7/7900HT (96-Well)	ViiA™ 7/7900HT (384-Well)	ViiA™ 7/7900HT Fast	310*	3100/3130 3700/3730
<b>96-Well Plates</b>																			
MicroAmp® Fast Optical 96-Well Reaction Plate (0.1 mL) - 10 plates	4346907	x																	
MicroAmp® Fast Optical 96-Well Reaction Plate with Barcode (0.1 mL) - 20 plates	4346906	x																	
MicroAmp® Fast Optical 96-Well Reaction Plate with Barcode (0.1 mL) - 200 plates	4366932	x																	
MicroAmp® Optical 96-Well Reaction Plate - 10 plates	N8010560		x			x	x	x				x	x		x				
MicroAmp® Optical 96-Well Reaction Plate - 500 plates	4316813		x			x	x	x				x	x		x				
MicroAmp® Optical 96-Well Reaction Plate with Barcode - 20 plates	4306737		x			x	x	x				x	x		x				
MicroAmp® Optical 96-Well Reaction Plate with Barcode - 500 plates	4326659		x			x	x	x				x	x		x				
MicroAmp® Optical 96-Well Reaction Plate with Barcode and Optical Caps - 20 plates	403012		x			x		x				x	x		x				
MicroAmp® Optical 96-Well Reaction Plate with Barcode and Optical Adhesive Films - 100 plates	4314320		x					x				x	x		x				
<b>384-Well Plates</b>																			
MicroAmp® Optical 384-Well Reaction Plate with Barcode - 50 plates	4309849			x												x			
MicroAmp® Optical 384-Well Reaction Plate with Barcode - 500 plates	4326270			x												x			
MicroAmp® Optical 384-Well Reaction Plate with Barcode - 1,000 plates	4343814			x												x			
MicroAmp® Optical 384-Well Reaction Plate - 1,000 plates	4343370			x												x			
<b>48-Well Plates</b>																			
MicroAmp® Fast Optical 48-Well Reaction Plate - 20 plates	4375816	x																	
<b>Single Tubes</b>																			
MicroAmp® Fast Reaction Tube with Cap (0.1 mL) - 1,000 tubes	4358297	x																	
MicroAmp® Reaction Tube with Cap (0.2 mL) - 1,000 tubes	N8010540		x			x		x											
MicroAmp® Reaction Tube with Cap (0.2 mL) - 10,000 tubes	N8011540		x			x		x											
MicroAmp® Reaction Tube with Cap, Assorted Colors (0.2 mL) - 1,000 tubes	N8010840		x			x		x											
MicroAmp® Reaction Tube with Cap, Autoclaved (0.2 mL) - 1,000 tubes	N8010612		x			x		x											
MicroAmp® Reaction Tube without Cap (0.2 mL) - 2,000 tubes	N8010533		x			x		x											
MicroAmp® Reaction Tube without Cap (0.2 mL) - 10,000 tubes	N8011533		x			x		x											
MicroAmp® Reaction Tube without Cap, Assorted Colors (0.2 mL) - 1,000 tubes	N8010833		x			x		x											
MicroAmp® Optical Tube without Cap (0.2 mL) - 2,000 tubes	N8010933		x			x		x					x	x					
GeneAmp® Thin-Walled Reaction Tube with Flat Cap (0.5 mL) - 1,000 tubes	N8010737				x														
GeneAmp® Thin-Walled Reaction Tube with Domed Cap (0.5 mL) - 2,000 tubes	N8010537				x														
GeneAmp® Thin-Walled Reaction Tube with Domed Cap, Autoclaved (0.5 mL) - 1,000 tubes	N8010611				x														
MicroAmp® 96-Well Reaction Tube/Tray/Retainer Set (0.2 mL) - 20 trays	403083					x		x											
MicroAmp® 96-Well Reaction Tube/Tray/Retainer Set (0.2 mL) - 100 trays	403086					x		x											
<b>Strip Tubes and Caps</b>																			
MicroAmp® 12-Cap Strip - 200 strips	N8010534	x	x			x		x											
MicroAmp® 12-Cap Strip - 1,000 strips	N8011534	x	x			x		x											
MicroAmp® 12-Cap Strip, Assorted Colors - 200 strips	N8010834	x	x			x		x											
MicroAmp® Fast 8-Tube Strip (0.1 mL) - 125 strips	4358293	x																	
MicroAmp® 8-Tube Strip (0.2 mL) - 125 strips	N8010580		x			x		x											
MicroAmp® 8-Tube Strip, Assorted Colors (0.2 mL) - 120 strips	N8010838		x			x		x											
MicroAmp® Optical 8-Tube Strip (0.2 mL) - 125 strips	4316567		x			x		x					x	x		x			
MicroAmp® 8-Cap Strip - 300 strips	N8010535	x	x			x		x											
MicroAmp® 8-Cap Strip - 1,500 strips	N8011535	x	x			x		x											
MicroAmp® 8-Cap Strip, Assorted Colors - 300 strips	N8010835	x	x			x		x											
MicroAmp® Optical 8-Cap Strip - 300 strips	4323032	x	x			x		x					x	x		x			
<b>Seals and Covers</b>																			
MicroAmp® Clear Adhesive Film - 100 films	4306311	x	x	x				x		x									
MicroAmp® Optical Adhesive Film - 25 films	4360954	x	x	x				x		x		x	x	x	x	x	x		
MicroAmp® Optical Adhesive Film - 100 films	4311971	x	x	x				x		x		x	x	x	x	x	x		
MicroAmp® 48-Well Optical Adhesive Film - 100 films	4375323	x																	
MicroAmp® 48-Well Optical Adhesive Film - 25 films	4375928	x																	
MicroAmp® Optical Adhesive Film Kit - 1 kit	4313663																		
MicroAmp® 96-Well Full Plate Cover - 5 covers	N8010550																		
<b>Reaction Trays</b>																			
MicroAmp® 96-Well Tray/Retainer Set - 10 trays	403081							x		x									
MicroAmp® 96-Well Tray - 10 trays	N8010541							x		x									
MicroAmp® 96-Well Tray for VeriFlex™ Blocks - 10 trays	4379983	x	x																
MicroAmp® 96-Well Tray/Retainer Set for Veriti® Systems - 10 trays	4381850		x																
MicroAmp® Fast 48-Well Tray - 10 trays	4375282																		
<b>Compression Pads</b>																			
MicroAmp® Optical Film Compression Pad - 5 pads	4312639																		
MicroAmp® Snap-On Optical Film Compression Pad - 9 pads	4333292																		

\* Requires the use of proper adapters.

For research use only. Not intended for any animal or human therapeutic or diagnostic use, unless otherwise stated.

NOTICE TO PURCHASER: PLEASE REFER TO THE USER'S GUIDE OR PACKAGE INSERT OF THE PRODUCTS NAMED HEREIN FOR LIMITED LABEL LICENSE OR DISCLAIMER INFORMATION.

© 2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. C016834 0111